



## quickAmp™ Plant Genomic DNA PCR Kit

### Storage

The kit should be stored at  $-20^{\circ}\text{C}$ .

### Procedure

All steps are performed at room temperature unless otherwise noted.

#### A. Plant Cell DNA Extraction

1. Rinse the paper punch and forceps in 70% ethanol prior to use and between the handling of different samples.
2. Punch a 0.5 to 0.7 cm disk of leaf tissue or frozen leaf tissue into a 1.5 ml or 2 ml microcentrifuge tube using a standard one-hole paper punch.
3. Pipette 100  $\mu\text{l}$  of Extraction Solution into the tube and vortex briefly.
4. Incubate samples at  $95^{\circ}\text{C}$  for 10 minutes.
5. Add 100  $\mu\text{l}$  of Neutralization Solution to sample and mix by vortexing.
6. Store the neutralized leaf extract at  $4^{\circ}\text{C}$  or use immediately in PCR. Continue with Section B, Step 1.

#### B. PCR Amplification

**ACTGene Taq mastermix, ACTGene Hot-Start Taq mastermix, ACTGene Blue/Red Taq mastermix, or ACTGene Blue/Red Hot-Start Taq mastermix** can be used for the following PCR amplification steps. Typical final primer concentrations are approximately 0.4  $\mu\text{M}$  each. Primer concentrations and thermal cycling parameters should be optimized depending on the system.

1. Add 10  $\mu\text{l}$  of PCR mastermix (2x), 4  $\mu\text{l}$  of leaf disk extract, and proper amount of primers to a thin-walled PCR microcentrifuge tube or plate. Bring the total volume to 20  $\mu\text{l}$  with PCR grade water. Gently mix well. Note: If less than 4  $\mu\text{l}$  of leaf disk extract is added to the PCR reaction volume, use a 50:50 mixture of Extraction Solution : Neutralization Solution to compensate the volume of tissue extract to 4  $\mu\text{l}$ .
2. Perform thermal cycling using optimized amplification parameters.
3. Load the amplified DNA onto agarose gel after the PCR is completed.