



quickAmp™ Tissue DNA Amplification Kit used to extract and amplify genomic DNA from human cells.

M DNA Ladder

- 1 Human hair follicle β -globin
- 2 Human buccal cell β -globin

Product Introduction

Product Description: quickAmp™ PCR kit includes reagents and enzymes for DNA genomic sample extraction and amplification, delivered in all in one convenient kit eliminating time consuming enzyme digestion, organic extraction and conventional centrifugation steps.

Typical Applications:
Genotyping, Transgene detection, Knockout analysis, and etc.

Product Highlights

- Quick Genomic PCR (<15 minutes) kit from genomic DNA extraction to amplification.
- Versatile kit for various sample types and applications:
 - Tissue (for mouse-tails, Hair, Animal Tissue, Saliva, & etc.)
 - Plant
 - Seed
 - Blood
- Sensitive only require small amount of starting material (see protocol for details by sample type).

Order Information

| | |
|---------------|---|
| A1500-100RXN | quickAmp™ Tissue DNA Amplification Kit, 100 rxns |
| A1500-1000RXN | quickAmp™ Tissue DNA Amplification Kit, 1000 rxns |
| A2500-100RXN | quickAmp™ Plant DNA Amplification Kit, 100 rxns |
| A2500-1000RXN | quickAmp™ Plant DNA Amplification Kit, 1000 rxns |
| A3500-100RXN | quickAmp™ Seed DNA Amplification Kit, 100 rxns |
| A3500-1000RXN | quickAmp™ Seed DNA Amplification Kit, 1000 rxns |
| A4500-100RXN | quickAmp™ Blood DNA Amplification Kit, 100 rxns |
| A4500-1000RXN | quickAmp™ Blood DNA Amplification Kit, 1000 rxns |

ACTGene, Inc.

P.O. Box 7011

Piscataway, NJ 08855

Tel: (973) 978-7888 | Fax: (732)469-0333

E-mail: sales@actgene.com